


NIMoSoft

Fetal RHD qPCR Analysis Software

The NIMoSoft software associated with the NIMoTest® Fetal RHD qPCR kit is the tool that will undoubtedly facilitate the analysis and interpretation of results.

It ensures complete traceability of actions. It automatically generates a configuration sheet and imports it into the qPCR device. It validates experiments and interprets results and integrates them into LIMS.



Standard ISO 15189

Complete Workflow

- Template creation 96 and 384 wells plate
- Template transfer to qPCR devices

Multiple qPCR devices validated

QuantStudio 7 - Life technologies
qTOWER3 - Analytic Jena AG
LC480 - Roche Life Science
M2000rt - Abbott



- Data analysis from qPCR devices
- Two level validation
- Export to LIMS

